

OXYGEN CONCENTRATOR

1L



SPECIFICATION

- Oxygen flow: 1-5 L/minute adjustable Flow
- Running Noise: ≤ 45 dB
- Dimension: 260x215x370(mm)
- Output Power: 90W AV(average)
- Caution: It shall not be used in surgery, emergency nor for critically ill patients
- Concentration: $93 \pm 3\%$
- Total Weight: 5.5KG
- Working Voltage: Alternating current 220V 50Hz ± 1 Hz
- Scope: Oxygen therapy and oxygen care for home and hospital.

