



DISPOSABLE MEDICAL PROTECTIVE CLOTHING

It is suitable for medical staff to contact with potentially infectious patients' blood, body fluid, particulate matter in the air to provide barrier and protection







PRODUCT PARAMETER



Product Name: Disposable medical protective clothing
Standard: GB19082-2009
Material: PE breathable film, non woven EVA rubber strip
Fabric: SS polypropylene non-woven fabric and PE breathable film 65gsm
Specification: One-piece(non sterile)160(S)/165(M)/170(L)/175(XL)/180(XXL)/185(XXXL)
Structure and composition: made of non-woven composite film, EVA rubber strip and zipper; designed with hooded tops and trousers.

PRODUCT DETAILS

	<p>Elastic Mouth Design</p> <p>The cuff trouser leg is designed with elastic opening for easy expansion.</p>
<p>Smooth Surface Without Moisture</p> <p>Can effectively prevent bacteria from adhering.</p>	
	<p>Sealing With Rubber Strip</p> <p>Effective isolation suture area zipper and other locations bacterial virus invasion.</p>
<p>Waist Elastic Design</p> <p>To meet the needs of different figures.</p>	

 <p>▪ Front</p>	 <p>▪ Side</p>
 <p>▪ Back</p>	

PACKAGING INFORMATION

