

10L Oxygen Concentrator

YS-501

An oxygen concentrator is a small oxygen generator would separate Oxygen from the air and concentrate to a high purity oxygen for inhalation.

YS-501 is a perfect product to be used at home or in hospitals. Many of our customers has been benefited from selling this item!



Features

- Reliable and durable oil free compressor;
- Large HD displayer, touch screen;
- 24 hours continuous working;
- Suitable for both home or hospitals;
- High oxygen purity 93%±3%;
- Dual flow for multiple users;

Specification

- ❖ Unit size: 335*320*620mm
- ❖ Power Supply: 220V/50Hz or 110V/60Hz
- ❖ Oxygen concentration: 93±3%;
- ❖ Noise Level: <55dB
- ❖ Flow Rate: 1- 10L/min flow rate optional
- ❖ Accuracy: Accurate to 0.1L/min
- ❖ Rating Current: 1.15A
- ❖ Dual Flow: Yes
- ❖ Nebulization: Yes
- ❖ Low Oxygen Concentration alarm: Yes
- ❖ Certificates: CE/ISO13485

Other functions

- ❖ Automatic alarm detection system
- ❖ Photosensitive large LED displayer
- ❖ Oxygen purity monitoring
- ❖ Cumulative timing
- ❖ Power failure alarm and other functions

Packing Detail

	Main unit	Oxygen tube	Nebulizer kit	Manual	Color box
Qty	1	1	1	1	1
1pcs/ctn N.W: 19.2kg G.W: 21.7kg Carton size: 535*340*680mm					

PROMED HEALTHCARE LIMITED

Add: Shenzhen | Hongkong | China

Email: SALES55@LLKMEDICAL.COM Web: www.llkmedical.com

PROMED 普尔迈