

# AIR COMPRESSION THERAPY SYSTEM

**PROMED 普尔迈**

IPC, with its multi-chamber sleeves squeezing patient legs repeatedly by air inflation and deflation, can substantially promote venous return, prevent deep vein thrombosis(DVT) and eliminate limb edema.



## FEATURE

- DVT&PTE Prevention
- Alleviate Limb Edema Rapidly
- Alleviate Lymphatic Swelling Rapidly
- Fat Dissolution
- Varicosity Improvement

## FEATURE

- Smart Venaflow detection to adjust therapy cycle
- Ergonomic design for easy operation
- 7" TFT touch screen user interface
- Auto sleeve detection
- Six grades adjustable inflation time
- 0-240 mmHg adjustable pressure for each chamber
- Replaceable lining design for sterilization
- Wire control for emergency stop
- Real time pressure monitoring with alarm protection
- 24 treatment modes including 18 customized ones



AIRPRO6000



AIPRO300



AIPRO600